

@551_EDB-DBCP

ANALYTE	METHOD	CONTAINER & PRESERVATIVES	HOLDING TIME	MRL	MDL	UNITS	LCS RANGE	MS RANGE	CAS #
Dibromochloropropane (DBCP)	EPA 551.1	60ml amber glass EDB1g(1%NaP/99%KP)+0.6%NH4CL	14 DAY	0.01	0.0084	ug/L	50 - 150	80 - 120	96-12-8
Ethylene Dibromide (EDB)	EPA 551.1	60ml amber glass EDB1g(1%NaP/99%KP)+0.6%NH4CL	14 DAY	0.01	0.0057	ug/L	50 - 150	80 - 120	106-93-4