

@ML531.2

ANALYTE	METHOD	CONTAINER & PRESERVATIVES	HOLDING TIME	MRL	MDL	UNITS	LCS RANGE	MS RANGE	CAS #
3-Hydroxycarbofuran	EPA 531.2	40ml amber glass vial 0.37g KH ₂ Citrate+6mg ThioSO ₄	28 DAY	0.5	0.143	ug/L	70 - 130	70 - 130	16655-82-6
Aldicarb (Temik)	EPA 531.2	40ml amber glass vial 0.37g KH ₂ Citrate+6mg ThioSO ₄	28 DAY	0.5	0.160	ug/L	70 - 130	70 - 130	116-06-3
Aldicarb sulfone	EPA 531.2	40ml amber glass vial 0.37g KH ₂ Citrate+6mg ThioSO ₄	28 DAY	0.5	0.172	ug/L	70 - 130	70 - 130	1646-88-4
Aldicarb sulfoxide	EPA 531.2	40ml amber glass vial 0.37g KH ₂ Citrate+6mg ThioSO ₄	28 DAY	0.5	0.142	ug/L	70 - 130	70 - 130	1646-87-3
Baygon	EPA 531.2	40ml amber glass vial 0.37g KH ₂ Citrate+6mg ThioSO ₄	28 DAY	0.5	0.247	ug/L	70 - 130	70 - 130	114-26-1
Carbaryl	EPA 531.2	40ml amber glass vial 0.37g KH ₂ Citrate+6mg ThioSO ₄	28 DAY	0.5	0.160	ug/L	70 - 130	70 - 130	63-25-2
Carbofuran (Furadan)	EPA 531.2	40ml amber glass vial 0.37g KH ₂ Citrate+6mg ThioSO ₄	28 DAY	0.5	0.104	ug/L	70 - 130	70 - 130	1563-66-2
Methiocarb	EPA 531.2	40ml amber glass vial 0.37g KH ₂ Citrate+6mg ThioSO ₄	28 DAY	0.5	0.114	ug/L	70 - 130	70 - 130	2032-65-7
Methomyl	EPA 531.2	40ml amber glass vial 0.37g KH ₂ Citrate+6mg ThioSO ₄	28 DAY	0.5	0.174	ug/L	70 - 130	70 - 130	16752-77-5
Oxamyl (Vydate)	EPA 531.2	40ml amber glass vial 0.37g KH ₂ Citrate+6mg ThioSO ₄	28 DAY	0.5	0.173	ug/L	70 - 130	70 - 130	23135-22-0