

@DIQUAT

ANALYTE	METHOD	CONTAINER & PRESERVATIVES	HOLDING TIME	MRL	MDL	UNITS	LCS RANGE	MS RANGE	CAS #
Diquat	EPA 549.2	1L amber poly no preservative	21 DAY	0.4	0.344	ug/L	70 - 130	70 - 130	85-00-7
Paraquat	EPA 549.2	1L amber poly no preservative	21 DAY	2	0.180	ug/L	70 - 130	70 - 130	1910-42-5