

@SPME

ANALYTE	METHOD	CONTAINER & PRESERVATIVES	HOLDING TIME	MRL	MDL	UNITS	LCS RANGE	MS RANGE	CAS #
Geosmin	SM 6040D	40ml amber glass vial SPME_80ul 32g/l NaOmadine	7 DAY	3	0.907	ng/L	75 - 125	70 - 130	16423-19-1
Methylisoborneol	SM 6040D	40ml amber glass vial SPME_80ul 32g/l NaOmadine	7 DAY	5	0.791	ng/L	75 - 125	70 - 130	2371-42-8