

## Glyphosate

<b>ANALYTE</b>	<b>METHOD</b>	<b>CONTAINER &amp; PRESERVATIVES</b>	<b>HOLDING TIME</b>	<b>MRL</b>	<b>MDL</b>	<b>UNITS</b>	<b>LCS RANGE</b>	<b>MS RANGE</b>	<b>CAS #</b>
Glyphosate	EPA 547	125ml amber glass no preservative	14 DAY	6	1.6	ug/L	70 - 130	70 - 130	1071-83-6