

@HAA5

ANALYTE	METHOD	CONTAINER & PRESERVATIVES	HOLDING TIME	MRL	MDL	UNITS	LCS RANGE	MS RANGE	CAS #
D/DBP Haloacetic Acids (HAA5)	SM 6251B	40ml amber glass vial 65mg NH4Cl	21 DAY	2	2	ug/L	80 - 120	70 - 130	-DBP (HAA5)
Dibromoacetic acid	SM 6251B	40ml amber glass vial 65mg NH4Cl	21 DAY	1	0.054	ug/L	80 - 120	84 - 122	631-64-1
Dichloroacetic acid	SM 6251B	40ml amber glass vial 65mg NH4Cl	21 DAY	1	0.105	ug/L	80 - 120	79 - 123	79-43-6
Monobromoacetic acid	SM 6251B	40ml amber glass vial 65mg NH4Cl	21 DAY	1	0.055	ug/L	80 - 120	81 - 122	79-08-3
Monochloroacetic acid	SM 6251B	40ml amber glass vial 65mg NH4Cl	21 DAY	2	0.410	ug/L	80 - 120	72 - 126	79-11-8
Trichloroacetic acid	SM 6251B	40ml amber glass vial 65mg NH4Cl	21 DAY	1	0.103	ug/L	80 - 120	82 - 124	76-03-9