

#PBCU

ANALYTE	METHOD	CONTAINER & PRESERVATIVES	HOLDING TIME	MRL	MDL	UNITS	LCS RANGE	MS RANGE	CAS #
Copper Total ICAP/MS	EPA 200.8	1L poly 4ml HNO3 (18%)	180 DAY	2	0.197	ug/L	85 - 115	70 - 130	7440-50-8
Lead Total ICAP/MS	EPA 200.8	1L poly 4ml HNO3 (18%)	180 DAY	0.5	0.038	ug/L	85 - 115	70 - 130	7439-92-1