

@ML505

ANALYTE	METHOD	CONTAINER & PRESERVATIVES	HOLDING TIME	MRL	MDL	UNITS	LCS RANGE	MS RANGE	CAS #
Alachlor (Alanex)	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.1	0.041	ug/L	70 - 130	65 - 135	15972-60-8
Aldrin	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.01	0.002	ug/L	70 - 130	65 - 135	309-00-2
Chlordane	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.1	0.032	ug/L	70 - 130	65 - 135	57-74-9
Dieldrin	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.01	0.005	ug/L	70 - 130	65 - 135	60-57-1
Endrin	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.01	0.005	ug/L	70 - 130	65 - 135	72-20-8
Heptachlor	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.01	0.003	ug/L	70 - 130	65 - 135	76-44-8
Heptachlor Epoxide	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.01	0.005	ug/L	70 - 130	65 - 135	1024-57-3
Lindane (gamma-BHC)	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.01	0.007	ug/L	70 - 130	65 - 135	58-89-9
Methoxychlor	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.05	0.022	ug/L	70 - 130	65 - 135	72-43-5
PCB 1016 Aroclor	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.08	0.022	ug/L	70 - 130	65 - 135	12674-11-2
PCB 1221 Aroclor	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.1	0.079	ug/L	70 - 130	65 - 135	11104-28-2
PCB 1232 Aroclor	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.1	0.085	ug/L	70 - 130	65 - 135	11141-16-5
PCB 1242 Aroclor	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.1	0.072	ug/L	70 - 130	65 - 135	53469-21-9
PCB 1248 Aroclor	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.1	0.023	ug/L	70 - 130	65 - 135	12672-29-6
PCB 1254 Aroclor	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.1	0.035	ug/L	70 - 130	65 - 135	11097-69-1
PCB 1260 Aroclor	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.1	0.033	ug/L	70 - 130	65 - 135	11096-82-5
Total PCBs	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.1	0.1	ug/L		65 - 135	1336-36-3
Toxaphene	EPA 505	40ml amber glass vial 1drop thio (8%)	8 DAY	0.5	0.083	ug/L	70 - 130	65 - 135	8001-35-2