

**@RA228 GA**

<b>ANALYTE</b>	<b>METHOD</b>	<b>CONTAINER &amp; PRESERVATIVES</b>	<b>HOLDING TIME</b>	<b>MRL</b>	<b>MDL</b>	<b>UNITS</b>	<b>LCS RANGE</b>	<b>MS RANGE</b>	<b>CAS #</b>
Radium 228	RA-228 GA	1L poly 4ml HNO3 (18%)	180 DAY	1		pCi/L	80 - 120	70 - 130	15262-20-1
Radium 228 Min Detect Activity	RA-228 GA	1L poly 4ml HNO3 (18%)	180 DAY			pCi/L			ium 228 MDA
Radium 228 Two Sigma Error	RA-228 GA	1L poly 4ml HNO3 (18%)	180 DAY			pCi/L			adium 228 Err