

#U-MST22

ANALYTE	METHOD	CONTAINER & PRESERVATIVES	HOLDING TIME	MRL	MDL	UNITS	LCS RANGE	MS RANGE	CAS #
Uranium ICAP/MS	EPA 200.8	500ml acid poly 2ml HNO3 (18%)	180 DAY	1	0.05698	ug/L	85 - 115	70 - 130	7440-61-1
Uranium by ICPMS as pCi/L	EPA 200.8	500ml acid poly 2ml HNO3 (18%)	180 DAY	0.7	0.022	pCi/L		70 - 130	7440-61-1092