

@ML515.4

ANALYTE	METHOD	CONTAINER & PRESERVATIVES	HOLDING TIME	MRL	MDL	UNITS	LCS RANGE	MS RANGE	CAS #
2,4,5-T	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.2	0.030	ug/L	70 - 130	70 - 130	93-76-5
2,4,5-TP (Silvex)	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.2	0.022	ug/L	70 - 130	70 - 130	93-72-1
2,4-D	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.1	0.028	ug/L	70 - 130	70 - 130	94-75-7
2,4-DB	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	2	0.524	ug/L	70 - 130	70 - 130	94-82-6
3,5-Dichlorobenzoic acid	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.5	0.069	ug/L	70 - 130	70 - 130	51-36-5
Acifluorfen	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.2	0.035	ug/L	70 - 130	70 - 130	50594-66-6
Bentazon	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.5	0.063	ug/L	70 - 130	70 - 130	25057-89-0
Dalapon	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	1	0.122	ug/L	70 - 130	70 - 130	75-99-0
Dicamba	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.1	0.017	ug/L	70 - 130	70 - 130	1918-00-9
Dichlorprop	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.5	0.056	ug/L	70 - 130	70 - 130	120-36-5
Dinoseb	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.2	0.024	ug/L	70 - 130	70 - 130	88-85-7
Pentachlorophenol	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.04	0.005	ug/L	70 - 130	70 - 130	87-86-5
Picloram	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.1	0.015	ug/L	70 - 130	70 - 130	01918-02-1
Tot DCPA Mono&Diacid Degradate	EPA 515.4	125ml amber glass 7mg SULFITE xls	14 DAY	0.1	0.044	ug/L	70 - 130	70 - 130	1861-32-1