

QTA[®] COTTON SOLUTIONS



Two minute analysis of cotton seed/meal without capital investment

Utilizing infrared technology, Eurofins QTA provides rapid, on-site testing solutions specifically for cotton seed and meal processors. QTA's expertise in model development delivers accurate and reliable results for multiple parameters - all in a 90 second test cycle.

Eurofins QTA is part of the Eurofins Scientific global network. With over 16,000 employees in more than 200 laboratories across 36 countries, Eurofins Scientific leads the world in the field of environmental laboratory testing and is a global market leader in agrosience, genomics, and central laboratory services.

Benefits of cotton analysis with Eurofins QTA

- **Fast and easy to use with minimal training**
- **Rapid analysis of multiple parameters**
- **Testing of Gases, Liquids, and Solids**
- **Nondestructive, no sample preparation**
- **Extensive reporting, CofA, and data services**
- **QTA handles all model validations and updates**
- **24/7 live technical support**

Material	Trait	Range
Cotton Meal	Moisture	10.5 – 15%
Cotton Meal	Oil	1 – 7.5%
Cotton Meal	Protein	30 – 46%
Cotton Seed	Total lint	1 – 11.6%
Cotton Seed	FFA	0.15-0.84%
Cotton Seed	Moisture	6.6-12.6%
Cotton Seed	Oil	17.2-23.5%
Cotton Hull	Oil	0.5 – 3.5%
Cotton Seed	Color	